

I'm not a bot



Clear blue ovu test

Trying for a baby? Double your chances of getting pregnant (1) with Clearblue Advanced Digital Ovulation kit. You can also track your cycles in our FREE app, designed for use with Clearblue Advanced Digital or Digital Ovulation Tests. Clearblue Advanced Digital Ovulation kit typically identifies 4 or more fertile days (2) when you have the best chances to get pregnant each cycle. By accurately tracking 2 key fertility hormones, estrogen and luteinising hormone, this ovulation test not only identifies 2 peak fertility days, but also additional high fertility days before that. This means you have more opportunities to get pregnant - at least twice as many as with any other ovulation test (2). Unlike basal body temperature (BBT) charting, ovulation predictor kits detect the rise in luteinising hormone (LH) that rises prior to ovulation to accurately detect the day BEFORE, and day of ovulation. Clearblue Advanced Digital Ovulation Test detects estrogen which rises prior to the LH surge to detect more fertile days. As sperm can survive in your body for up to 5 days, you can become pregnant by having sex during the days leading up to ovulation (on high and peak fertility days) to maximise your chances of getting pregnant. The clear digital display shows a Clearblue smiley face to identify your best days to get pregnant. Clearblue Advanced Digital Ovulation Test is over 99% accurate at detecting the LH surge. (1) In first cycle of use, vs not using ovulation tests. (2) In a study of 87 women, 4 or more fertile days were identified in 80% of cycles using actual cycle length (2012). To find out if your phone is compatible and if the Clearblue Period & Cycle Tracker app is available in the country you live in, check out: clearblue.com/connectivity. Double your chances of getting pregnant in the first cycle of use, vs not using ovulation tests. Clear results, unique Clearblue 'smiley face' appears on your fertile days: Have sex on High and Peak Fertility days to maximise your chances of getting pregnant. Typically identifies 4 or more fertile days (1) each cycle. It tracks 2 key fertility hormones to identify your High AND Peak fertility days, so you can plan ahead and have more opportunities to get pregnant. No. 1 UK Doctor recommended brand in a survey of doctors in the UK (data on file). Providing pregnancy tests, ovulation kits and fertility tests for women for over 30 years. Clearblue Advanced Digital Ovulation Tests adapts to your unique hormone profile for cycle understanding. Over 99% accurate at detecting the LH surge. Easy to use vs ovulation test strips - 7 in 10 [73%] women did not find strip tests easy to use. (2) In a study of 111 women in the UK using and reading different test types. Data on file. Try FREE Period & Cycle Tracker app: view where you are in your cycle and see how the probability of peak fertility varies for women with the same cycle length as you. Use Clearblue Ovulation Tests to identify your peak days. Always read the manufacturer's instructions for any medication you are taking before testing. You may get misleading results if you are taking fertility drugs containing luteinising hormone (LH) or human Chorionic Gonadotrophin, or are taking antibiotic containing tetracyclines. Some fertility treatments such as clomiphene citrate may give misleading High Fertility results but Peak Fertility results should not be affected. If you have recently stopped using hormonal contraception your cycles may be irregular so you may wish to wait until you have had 2 cycles before testing. As prenatal care is very important for a baby's health, we recommend that you consult your doctor before you try to conceive. Check with your doctor if you are taking any medication or have any medical condition before planning a pregnancy. If you have a medically diagnosed fertility problem, ask your doctor if this test is suitable for you and if you do get unexpected results discuss them with your doctor. If you are under 35 years and haven't become pregnant, we recommend you see your doctor after 12 months. If you're over 35 years see your doctor if you haven't become pregnant after 6 months, and straight away if you are over 40. This IVD digital device meets the missions and immunity requirements of EN 61326-2-6. The EMC countermeasures employed within the electronic instrument will provide reasonable protection against electromagnetic interference effects likely to be encountered in the home environment. The following preventive warnings apply to EN 61326-2-6 compliant equipment. a) Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results. b) Do not use this instrument in close proximity to sources of strong electromagnetic radiation (e.g., mobile phones), as these may interfere with the proper operation. SPD Swiss Precision Diagnostics GmbH, Route de St Georges 47, 1213 Petit-Lancy, Geneva, Switzerland. GB: 0800 917 2710; IE: 1800 812 607. www.clearblue.com. Name & Address of Authorised Representative: MDSS GmbH, Schiffgraben 41 Hannover 30175, Germany; info@mdssar.com. Country of Origin: Switzerland. Address: Manufacturer Address SPD Swiss Precision Diagnostics GmbH, Route de St Georges 47, 1213 Petit-Lancy, Geneva, Switzerland. Consumer Care Return Address Customer Support, SPD Development Company Ltd, Clearblue Innovation Centre, Priory Business Park, Bedford, MK44 3UP United Kingdom